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Harlequin® UTI Chromogenic Agar (NCM1013)

Intended Use
Urinary Tract Infection (UTI) Chromogenic Agar is a chromogenic medium for the presumptive identification and confirmation of microorganisms causing urinary tract infections.

Description

UTI Chromogenic Agar contains a peptone mixture to provide nitrogen, vitamins, minerals and amino acids essential for growth. The medium includes two chromogenic substrates which are hydrolyzed by enzymes produced by *Enterococcus* spp, *Escherichia coli* and coliforms. It also includes phenylalanine and tryptophan providing a presumptive indication of the tryptophan deaminase activity, indicating the presence of Proteus spp., *Morganella* spp, and *Providencia* spp. (brown colonies).

One chromogen is metabolized by β-glucosidase enzyme activity, allowing the specific detection of enterococci which form blue or turquoise colonies. The other chromogen is hydrolyzed by β-galactosidase, an enzyme produced by *E. coli* which grows as pink colonies (in case of unreliable colony results, carry out Indol test). When bacteria utilize both chromogenic substrates, this results in dark blue - purple colonies, characteristic of coliform bacteria such as *E. aerogenes*, *K.pneumoniae* and *C. freundii*.

Typical Formulation

Peptone Mixture	16.0 g/L
Growth Factors	13.0 g/L
Tryptophan	2.0 g/L
Chromogenic Substrate	0.50 g/L
Agar	16.0 g/L

Final pH: 7.2 ± 0.2 at 25°C
Formula may be adjusted and/or supplemented as required to meet performance specifications.

Precaution
1. Refer to SDS

Preparation
1. Suspend 47.5 grams of the medium in one liter of distilled water.
2. Heat with frequent agitation and boil for one minute to completely dissolve the medium.
3. Autoclave at 121ºC for 15 minutes.
4. Cool to 45-50ºC.

Test Procedure
Consult appropriate references for the isolation and identification of anaerobic bacteria.

Quality Control Specifications
Dehydrated Appearance: Powder is homogeneous, free flowing and beige.

Prepared Appearance: Prepared medium is amber in color and slightly opalescent.

Minimum QC:

Escherichia coli ATCC 25922	Good growth	Pink colonies
Enterobacter aerogenes ATCC 13048	Good growth	Dark Blue colonies
Proteus miriabilis ATCC 13315	Good growth	Light Brown colonies
Staphylococcus aureus ATCC 25923	Good growth	White-Cream colonies

Effective Date: 1/17/2019
Revision: 0

620 Lesher Place • Lansing, MI 48912
800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com

Salmonella typhimurium ATCC 14028 Good growth Amber colonies

Results

Organism	Growth	Colony Color
Escherichia coli	Good	Pink
Enterobacter aerogenes	Good	Dark blue
Klebsiella pneumonieae	Good	Dark blue
Proteus miriabilis	Good	Light brown
Staphylococcus aureus	Good	White/cream (natural pigmentation)
Enterococcus faecalis	Good	Light blue
Pseudomonas aeruginosa	Good	Amber
Salmonella typhi	Good	Amber
Salmonella typhimurium	Good	Amber

Expiration

Refer to expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing or appearance has changed from the original color. Expiry applies to medium in its intact container when stored as directed.

Limitations of the Procedures

Chromogenic agars for UTI can give atypical results when atypical strains are present in the sample. For example, 45% of *Enterobacter cloacae* do not contain β-glucosidase, therefore resulting in pink colonies not distinguishable from *E. coli*. For confirmation, the Indol test must be performed.

Storage

Store dehydrated culture media at 2-25°C away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

References

Samra Z, Heifetz M, Talmor J, Bain E and Bahar J. Evaluation of Use of a New Chromogenic Agar in Detection of Urinary Tract Pathogens. Journal of Clinical Microbiolology. 1998;36 (4): 990-4.

Effective Date: 1/17/2019

Revision: 0

620 Lesher Place • Lansing, MI 48912
800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com