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Product 9822

Read the instructions carefully before starting the test.

NEOGEN®

ANSR®

for *E. coli* O157:H7

PERFORMANCE TESTED
AOAC RESEARCH INSTITUTE
LICENSE NUMBER 111502

INTENDED USE

In an AOAC Research Institute Performance Tested℠ Method Study, ANSR® for *E. coli* O157:H7 was found to be an effective procedure for the detection of *E. coli* O157:H7 in raw ground beef, raw beef trim, spinach, and sprout irrigation water.

ASSAY PRINCIPLES

ANSR for *E. coli* O157:H7 is an isothermal, amplified nucleic acid assay. The ANSR for *E. coli* O157:H7 method is based on nicking enzyme amplification reaction (NEAR™) technology. DNA is amplified through a mechanism of polymerization from the ends of nicks created in double-stranded DNA by the action of a specific endonuclease. Amplified target sequences are detected in real time using fluorescent molecular beacon probes.

A two-stage lysis reaction is performed, first at 37 ± 2°C for ten minutes, then at 80 ± 2°C for 20 minutes. Next, a portion of the lysed sample is transferred to a strip tube containing lyophilized ANSR reagents. The tubes are sealed and incubated at 56 ± 2°C on the ANSR reader. Results are generated by the reader and displayed in the ANSR software within ten minutes. Positive results may be confirmed from the enrichment cultures following standard procedures. Each tube of ANSR reagents contains an internal positive control, ensuring that the reagents are functioning properly.

INTENDED USER

The ANSR® for *E. coli* O157:H7 test is designed for use by personnel with appropriate training in microbiology. Training in the use of the ANSR test system is available through NEOGEN®.

MATERIALS PROVIDED

1. 12 strips of 8 cluster tubes, 1.2 mL
2. 12 strips of 8 reaction tubes, 200 µL, containing lyophilized ANSR for *E. coli* O157:H7 reagents in 2 sealed foil pouches with desiccant pack
3. 12 strips of 8 permanent caps for reaction tubes
4. 1 bottle lysis reagent suspension buffer, 60 mL
5. 3 vials containing lyophilized lysis reagents
6. 1 kit insert

EQUIPMENT REQUIRED

1. ANSR reader (NEOGEN item 9828)
2. Computer and software for connection to ANSR reader (NEOGEN item 9832)
3. 1 dual heater block with aluminum block inserts for 1.2 mL cluster tubes, 80 ± 2°C and 37 ± 2°C (NEOGEN items 9386-48D, 9829-48) or 2 single heater blocks with aluminum block inserts for 1.2 mL cluster tubes, 80 ± 2°C and 37 ± 2°C (NEOGEN items 9386-48D, 9829-48)
4. Pipettor, 100–1000 µL (NEOGEN item 9463)
5. Pipette tip rack, 100–1000 µL, sterile, 96 tips (NEOGEN item 9487)
6. Pipettor, 10–100 µL, 8-channel (NEOGEN item 9388)
7. Pipette tips, 100 µL, sterile, filtered, 96 tips (NEOGEN item 9389)
8. Vortex mixer with an adjustable speed (NEOGEN item 9494)
9. Stomacher (optional)
10. 3 thermometers, traceable (NEOGEN item 9518)
11. Timer, 3-channel (NEOGEN item 9426)
12. Optional-for-use heater block with 0.2 mL reaction tube aluminum block insert, 56 ± 2°C (NEOGEN items 9386-48D, 9829-48)
13. Webcam (NEOGEN item WEBCAM)
14. ANSR ethernet cable (NEOGEN item 9835)
15. 10 mL pipette pump (NEOGEN item 9277)
16. Pipettes, sterile serological 10 mL (NEOGEN item 8686)
17. 40-slot, 20 mm test tube rack, autoclavable (NEOGEN item 9553)
18. Pipettor, 20–200 mL (NEOGEN item 9276)

OTHER MATERIALS REQUIRED

1. Stomacher-type bags for sample enrichment. Filtered bags are recommended (NEOGEN item 6833)
2. Graduated cylinder, 250 mL (NEOGEN item 9368)

MEDIA ENRICHMENT BROTH REQUIRED

1. ANSR for *E. coli* Enrichment Broth (NEOGEN item 9815)
2. Tryptic Soy Broth, Modified with Acid Digest of Casein (mTSB) (NEOGEN item 7167A)
3. 1L of purified water
4. Phosphate Buffered Saline (NEOGEN item 8425)

STORAGE

Store ANSR® reagents at 2–8°C. After removing the reaction tubes from the foil pouch, promptly re-seal the pouch. Leave the desiccant pack in the pouch at all times.

PRECAUTIONS

1. Use good microbiology laboratory practices.
2. Dispose of used pipette tips in a covered container containing a fresh solution of 10% bleach. The 10% bleach solution should be made fresh each day. Use 9 parts water with 1 part commercially available household-strength bleach to make the 10% bleach solution. Undiluted stock solutions of bleach should be used within 30 days after opening.
3. Discard bleach solution and tips as regular waste at the end of each day.
4. Do not use reagents beyond the expiration date.
5. Enrichment media should be prewarmed to incubation temperature (all protocols) before use.
6. Use of enrichment media and incubation times or temperatures other than those specified may lead to inaccurate results.
7. Remove the reaction tubes from the foil pouch just before use and keep covered until the heating process begins. Reseal the pouch containing the remaining reaction tubes to avoid prolonged exposure to light. More than 15 minutes of total exposure time may lead to incorrect results.
8. Do not, under any circumstance, remove the caps from reaction tubes after the assay has been started. This is essential in order to prevent accidental contamination of the environment with amplification products.
9. Exercise care in all pipetting steps to avoid cross-contamination of samples.
10. Complete all assay steps in sequence, avoiding delays between steps.
11. Tap the reaction tubes on bench top to make sure the lyophilized reagents are at the bottom of the tube prior to adding the lysed sample.

PREPARATION OF ENRICHMENT BROTH

ANSR for *E. coli* Enrichment Broth

Rehydrate 39.6 g of the medium in 1 L of sterile water prewarmed to 42 ± 1°C. Use rehydrated media the same day as prepared. Do not autoclave.

mTSB

Dissolve 43 g of the medium into 1 L of purified water. Autoclave at 121 ± 1°C for 15 minutes.

SAMPLE ENRICHMENT

Sprout irrigation water

1. Weigh out 125 mL of sprout rinse water, in a Stomacher-type bag.
2. Add 1.125 L ANSR for *E. coli* Enrichment Broth prewarmed to 42 ± 1°C to the bag. Swirl to mix.
3. Incubate the culture at 42 ± 1°C for 12–24 hours.

Leafy green vegetables

1. Weigh out 200 g of sample in a Stomacher-type bag.
2. Add 600 mL ANSR for *E. coli* Enrichment Broth prewarmed to 42 ± 1°C to the bag. Swirl to mix.
3. Incubate the culture at 42 ± 1°C for 24–26 hours.

Raw ground beef

1. Weigh out 325 g of sample in a Stomacher-type bag.
2. Add 1.625 L of ANSR® for *E. coli* Enrichment Broth or 975 mL mTSB prewarmed to 42 ± 1°C to the bag. Pummel, blend or hand massage until clumps are dispersed.
3. Incubate the culture at 42 ± 1°C for 12–24 hours for ANSR *E. coli* Enrichment Broth or 18–20 hours for mTSB.

Raw beef trim

1. Weigh out 375 g of sample in a Stomacher-type bag.
2. Add 1.875 L of ANSR for *E. coli* Enrichment Broth prewarmed to 42 ± 1°C to the bag. Pummel, blend or hand massage until clumps are dispersed.
3. Incubate the culture at 42 ± 1°C for 12–24 hours.

Raw beef trim (USDA Reference Method Enrichment)

1. Weigh out 25 g of sample in a Stomacher-type bag.
2. Add 225 mL of mTSB prewarmed to 42 ± 1°C to the bag. Pummel, blend, or hand massage until clumps are dispersed.
3. Incubate the culture at 42 ± 1°C for 18–20 hours.
 Note: When using smaller sample sizes than those described above, maintain the same ratio of enrichment volume to sample weight and incubation conditions.

LYSIS REAGENT SOLUTION PREPARATION

Reconstitute 1 vial of the lyophilized lysis reagents with 18 mL of the lysis reagent suspension buffer by adding the buffer to the reagent vial. Swirl gently to mix.
Note: 1 vial of the lysis reagents is enough for approximately 32 samples. Prepared lysis reagent solution can be stored at 2–8°C for up to 30 days.

ANSR TEST PROCEDURE

Prior to starting the assay:

1. Preheat 1 lysis heater block to 80 ± 2°C. Preheat the second lysis heater block to 37 ± 2°C. If using the optional single heater, preheat to 56 ± 2°C. Use the thermometer for the temperature reading.
2. Remove the foil pouch containing the reaction tubes from the refrigerator and allow them to warm at room temperature for 15 minutes. To avoid excess light exposure, leave the reaction tubes in the foil pouch until they are needed.
 Note: The lysis reagent solution should be at room temperature for 15 minutes before use.
3. Connect the ANSR reader to the computer via USB or ethernet and turn the computer on.
4. Turn on the ANSR reader. The reader will preheat to 56 ± 2°C.
5. Start the ANSR software and click the connect button. Input the sample IDs, lot number, and user information.
 Note: For instructions on using the reader and software, see the user guide that came with the ANSR reader.

ASSAY PROCEDURE

Note: Enrichment cultures that are highly turbid or highly colored may require a 1:10 dilution in PBS before lysing the sample. Contact NEOGEN® Technical Services for more information. The AOAC approved matrices do not require dilution.

OPTIONAL: Pre-pierced plate sealing film (BIO RAD Item 3600040) can be placed over the 1.2 mL cluster tubes if desired.

1. Add 50 µL of the enrichment culture to each 1.2 mL cluster tube using a pipettor with 100 µL filtered tips.
 Note: Cluster tubes may be pulled apart to provide the number of tubes needed.
2. Add 450 µL of the lysis reagent solution to each cluster tube containing culture. Be sure to switch pipette tips between samples.
 Note: Return the lysis reagent solution to the refrigerator within 1 hour of its use.
3. Incubate each cluster tube at 37 ± 2°C for 10 minutes.
4. Immediately transfer each cluster tube to the 80 ± 2°C heater block and incubate for 20 minutes.
 Note: The 80 ± 2°C incubation time may be extended to a total of 1 hour for the purpose of managing staggered assay start times.
5. For 3 to 5 minutes before the end of the lysis step, preheat the ANSR® reagents to 56 ± 2°C by placing each reaction tube in the ANSR reader.
 Optional: A separate heat block can be used. It should be heated to 56 ± 2°C.
 Note: The strip of reaction tubes may be cut to provide the number of tubes needed. Keep all unused tubes in the sealed foil pack. Ensure the pellet in the reaction tube is at the bottom by tapping the tubes gently on the bench top.
6. After the completion of the 20–60 minute lysis incubation, remove and discard the caps from each reaction tube in the ANSR reader.
 Important: Proceed with steps 7–9 without delay. The transfer of the sample from the lysis tubes to the reaction tubes should be completed within 1 minute.
7. Using an 8-channel pipette and 100 µL filter tips, carefully transfer 50 µL from the top third of each lysed sample in the cluster tube to each reaction tube. The assay's performance will be interfered if there's debris accumulation at the bottom of the lysis tube(s). Avoid the transfer of debris by aspirating from the top third of lysis tube. Do not prime the pipette tips and do not mix before aspirating.
 Place the provided permanent caps on the reaction tubes.
 Note: Lysed sample may be transferred from the same cluster tube a maximum of 3 times.
8. Remove the strips(s) of tubes from the reader (or 56 ± 2°C heat block if 1 was used) and vortex briefly (about 2 seconds), at medium or high speed (2500 RPM). Quickly visually check each reaction tube to assure that no bubbles are on the bottom or in the middle of the sample. A quick tap of the tubes should release any bubbles from the bottom or middle to the top. Then place into the reader without delay. Close the reader's lid.

Note: The reader will not provide accurate results if the lid is open. Keep the lid closed at all times while the assay is running. Contamination may occur if the permanent caps are not placed on the reaction tubes or if the permanent caps are removed.

9. Click start in the ANSR® software to begin the 10 minute assay.
10. Results will be displayed as positive, negative, or invalid once the assay is finished. If the result displays an invalid, the test must be repeated. A single repeat test can be run starting at step 6. Alternatively, start over from step 1 if the lysis from step 6 has been at 80 ± 2°C for more than 60 minutes.

INTERPRETATION OF RESULTS

The ANSR software will indicate the test results as positive or negative for the presence of *E. coli* O157:H7 in the enriched sample. In addition, the real-time fluorescence curve generated from the assay can be viewed.

CONFIRMATION

Samples producing positive ANSR results may be confirmed by streaking the enrichment cultures to *E. coli* O157:H7 selective/differential agar media and continuing with identification of presumptive colonies following standard methods [1,2]. Immunomagnetic separation using magnetic particles coated with anti-O157:H7 antibodies is recommended before to plating. CHROMagar O157 agar is recommended as the plating medium.

DISPOSAL

Enrichment cultures and used lysis tubes should be disposed of as biohazard waste. The preferred method of treatment for biohazard waste is autoclaving. Items that cannot be autoclaved, such as used reaction tubes and caps, should be immersed in a fresh 10% bleach solution that is made daily. Consult with the safety advisor of your facility for detailed instructions.

Do not remove the permanent caps, for any reason, from the ANSR reaction tubes once the assay has started — even when disposing of them. The reaction tubes can be disposed of as non-biohazardous waste. It is recommended that they be placed in a sealable container or plastic bag with 10% bleach and immediately disposed of to protect against accidental opening.

CUSTOMER SERVICE

NEOGEN® Customer and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on our website at NEOGEN.com, or by calling us at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For NEOGEN®'s full terms and conditions, please visit NEOGEN.com/terms-and-conditions.

WARRANTY

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

NEOGEN ANSR Molecular Diagnostic for Foodborne Pathogen Detection – limited use label license

SYTO® 82

SYTO® 82 contained within this product is provided under an intellectual property license from Life Technologies Corporation, Eugene, OR, and may be used for *in vitro* detection and analysis of (i) food, feeds and beverages, including nutraceuticals, (ii) ingredients for food, feeds and beverages, (iii) process samples from food, feed and beverage preparation, distribution and delivery, and (iv) water from any source for human consumption, all for the purpose of safety and quality assurance. The buyer must not sell or otherwise transfer this product or its components for any other use, including but not limited to: human *in vitro*, veterinary, identity or paternity testing, forensics, or *in vivo* detection of nucleic acid sequences in living beings, or cells. For information on purchasing a license for SYTO 82 for purposes other than food, beverage and water safety, and quality assurance, contact Life Technologies Corporation at outlicensing@lifetech.com.

Molecular beacon probes

One or more molecular beacon probes contained within this product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for tests on food products, feeds, beverages, and water.

NEAR™ Technology

This product utilizes the patented NEAR™ isothermal technology and is sold under license from Ionian Technologies, San Diego, CA, and may be used under Ionian Technologies patent rights only for tests on food, beverage, and water safety.

References

[1] USDA-FSIS (2015) Microbiology Laboratory Guidebook, chapter 5.09
http://www.fsis.usda.gov/wps/wcm/connect/51507fdb-dded-47f7-862d-ad80c3ee1738/MLG-5.pdf?MOD=AJPERES
[2] USDA-FDA (2014) Bacteriological Analytical Manual, chapter 4A
http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm070080.htm

PRODUCTS AVAILABLE

9465 **Environmental sampling sponges** (with Dey Engley (DE) neutralizing broth) — pack of 100

36002 **Sampling sponges, pre-moistened, HiCap** — pack of 100

9822 **ANSR for *E. coli* O157:H7** — up to 96 tests

9824 **ANSR for *Listeria monocytogenes*** — up to 96 tests

9870 **ANSR for *Salmonella*** — up to 96 tests

9871 **ANSR for *Listeria*** — up to 96 tests

9872 **ANSR for *Campylobacter*** — up to 96 tests

9873 **ANSR for *Listeria* Right Now™** — up to 96 tests

9837 **Complete ANSR Equipment Kit** — reader, computer, pipette set, dry bath heaters and inserts, thermometers, vortex, time, rack

North America
NEOGEN Headquarters
800.234.5333 (U.S./Canada)
foodsafety@neogen.com
NEOGEN.com

Europe, Middle East, and Africa
NEOGEN Europe
+ 44 (0) 1292 525 600
info_uk@neogeneurope.com
NEOGEN.com

Mexico
NEOGEN Latinoamerica
+52 (55) 5254.8235
informacion@neogenlac.com
NEOGEN.com

Brazil
NEOGEN do Brasil
+55 19 3935.3727
info@neogendobrasil.com.br
NEOGEN.com

China
NEOGEN Bio-Scientific Technology
+86 21 6271 7013
info@neogenchina.com.cn
www.neogenchina.com.cn

India
NEOGEN Food and Animal Security
+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

©NEOGEN Corporation, 2020. NEOGEN and ANSR are registered trademarks of NEOGEN Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.
Patent: https://www.neogen.com/patents/

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ANSR®

para *E. coli* O157: H7

PERFORMANCE TESTED AOAC RESEARCH INSTITUTE
LICENSE NUMBER 111502

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