SNAPSHOOT
ENGINEERED FOR A SAFER TOMORROW

현장을 이해하고, 필요한 솔루션을 연결합니다.

산업안전 · 추락방지 시스템
산업용 소재 · 건축자재
제품 제안부터 공급·시공·설치·교육까지



COMPANY INFORMATION

스냅슛 (SNAPSHOOT)

ADDRESS (15431) 경기도 안산시 단원구 풍전로 37-9 트리플렉스 301동 118호
TEL 0507-1374-9229
FAX 0504-481-6749
E-MAIL bkkim@snapshoot.co.kr

TECHNICAL DOCUMENTS

인증서 · 데이터시트 · 기술자료

제품의 규격과 성능, 인증 및 적용 조건을 확인하기 위한 자료입니다.

ABOUT SNAPSHOOT

제품 공급을 넘어, 현장에 맞는 해결책으로

스냅슛은 산업현장의 작업 환경과 요구 조건을 바탕으로
적합한 제품과 시스템을 제안합니다. 현장 실사를 통해 설치 조건을 확인하고,
공급부터 시공·설치 및 사용 교육까지 필요한 과정을 함께합니다.

추락방지 시스템

안전대부착설비 · 수평·수직 라이프라인
레일 시스템 · 리지드 포스트

밀폐공간 진입·구조 솔루션

데빗암 · 삼각대 · 구조용 윈치
윈치 안전블록 · 구조장비

산업안전 보호구

전신형 안전대 · 랜야드 · 안전블록
작업 환경에 맞춘 개인보호구

산업용 소재·건축자재

산업용 테이프 · 접착제 · 기능성 소재
구조용 스크류 · 목조건축 자재

SNAPSHOOT
ENGINEERED FOR A SAFER TOMORROW

홈페이지 바로가기 →
문서 라이브러리 바로가기 →
시공사례 바로가기 →
문의하기 →

본 페이지는 스냅슛의 자료 제공 안내용 표지입니다. 첨부된 인증서 및 기술자료의 발행 주체, 인증 대상과 적용 범위는 해당 원문을 기준으로 확인해 주세요.

Customer Service
NEOGEN Customer and Technical Services can be contacted through NEOGEN.com and product training is available by request.

SDS Information
No SDS is required or available for this kit.

Terms and conditions
NEOGEN's full terms and conditions are available online.

Warranty
NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

NEOGEN.com

©NEOGEN Corporation, 2021. NEOGEN, and AccuPoint are registered trademarks of NEOGEN Corporation, Lansing, MI 48912 U.S. All other brand and product names are trademarks or registered trademarks of their respective companies. Patent: https://www.NEOGEN.com/patents/

Product Number: 9912

AccuPoint® Advanced Next Gen

Electronic Control

Intended Use
The AccuPoint® Advanced Next Generation (NG) Electronic Control Kit consists of two controls (positive and negative) designed to verify the calibration of AccuPoint Advanced NG instrument. The positive control is designed to emit a low level LED light, traced to a NIST reference material. The negative control is used to check for possible background light, and should yield a 0 reading.

Procedure
1. Record the serial number and acceptable relative light units (RLU) range for the positive control in the table below.
2. Power on the AccuPoint Advanced NG instrument.
3. Open the chamber door. Insert the negative control and close the door.

NEOGEN® NEOGEN.com

Customer Service
NEOGEN Customer and Technical Services can be contacted through NEOGEN.com and product training is available by request.

SDS Information
No SDS is required or available for this kit.

Terms and conditions
NEOGEN's full terms and conditions are available online.

Warranty
NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

NEOGEN.com

©NEOGEN Corporation, 2021. NEOGEN, and AccuPoint are registered trademarks of NEOGEN Corporation, Lansing, MI 48912 U.S. All other brand and product names are trademarks or registered trademarks of their respective companies. Patent: https://www.NEOGEN.com/patents/

17751A Official_9912_AccuPoint-Advanced-NG-Elec-Control_Kit-Insert

Product Number: 9912

AccuPoint® Advanced Next Gen

Electronic Control

Intended Use
The AccuPoint® Advanced Next Generation (NG) Electronic Control Kit consists of two controls (positive and negative) designed to verify the calibration of AccuPoint Advanced NG instrument. The positive control is designed to emit a low level LED light, traced to a NIST reference material. The negative control is used to check for possible background light, and should yield a 0 reading.

Procedure
1. Record the serial number and acceptable relative light units (RLU) range for the positive control in the table below.
2. Power on the AccuPoint Advanced NG instrument.
3. Open the chamber door. Insert the negative control and close the door.

NEOGEN® NEOGEN.com

4. Record the result in the table below.
5. Repeat for a total of 10 tests.
6. Repeat the process for the positive control to collect a total of 10 readings.
7. If the %CV is greater than 10%, or the readings are outside the acceptable RLU range, the reader is out of calibration. Contact Technical Services for a replacement.

	Positive Control	Negative Control
Serial Number		
Acceptable RLU Range		0–5 RLU
Test 1		
Test 2		
Test 3		
Test 4		
Test 5		
Test 6		
Test 7		
Test 8		
Test 9		
Test 10		
Average		
SD		
%CV (SD + average)		

Shelf Life

The positive and negative control samplers have capacity to perform 5,000 readings or a shelf life of 3 years from manufactured date.

Storage

When not in use, store the two controls in the protective case.

Interpretation of Results

Calculate the average, standard deviation and %CV (SD ÷ average) for the positive control readings. If results are outside the acceptable RLU range, or if the %CV is >10% for the positive control, contact your NEOGEN representative.

4. Record the result in the table below.
5. Repeat for a total of 10 tests.
6. Repeat the process for the positive control to collect a total of 10 readings.
7. If the %CV is greater than 10%, or the readings are outside the acceptable RLU range, the reader is out of calibration. Contact Technical Services for a replacement.

	Positive Control	Negative Control
Serial Number		
Acceptable RLU Range		0–5 RLU
Test 1		
Test 2		
Test 3		
Test 4		
Test 5		
Test 6		
Test 7		
Test 8		
Test 9		
Test 10		
Average		
SD		
%CV (SD + average)		

2

Shelf Life

The positive and negative control samplers have capacity to perform 5,000 readings or a shelf life of 3 years from manufactured date.

Storage

When not in use, store the two controls in the protective case.

Interpretation of Results

Calculate the average, standard deviation and %CV (SD ÷ average) for the positive control readings. If results are outside the acceptable RLU range, or if the %CV is >10% for the positive control, contact your NEOGEN representative.

3