SNAPSHOOT
ENGINEERED FOR A SAFER TOMORROW

현장을 이해하고, 필요한 솔루션을 연결합니다.

산업안전 · 추락방지 시스템
산업용 소재 · 건축자재
제품 제안부터 공급·시공·설치·교육까지



COMPANY INFORMATION

스냅슛 (SNAPSHOOT)

ADDRESS (15431) 경기도 안산시 단원구 풍전로 37-9 트리플렉스 301동 118호
TEL 0507-1374-9229
FAX 0504-481-6749
E-MAIL bkkim@snapshoot.co.kr

TECHNICAL DOCUMENTS

인증서 · 데이터시트 · 기술자료

제품의 규격과 성능, 인증 및 적용 조건을 확인하기 위한 자료입니다.

ABOUT SNAPSHOOT

제품 공급을 넘어, 현장에 맞는 해결책으로

스냅슛은 산업현장의 작업 환경과 요구 조건을 바탕으로
적합한 제품과 시스템을 제안합니다. 현장 실사를 통해 설치 조건을 확인하고,
공급부터 시공·설치 및 사용 교육까지 필요한 과정을 함께합니다.

추락방지 시스템

안전대부착설비 · 수평·수직 라이프라인
레일 시스템 · 리지드 포스트

밀폐공간 진입·구조 솔루션

데빗암 · 삼각대 · 구조용 윈치
윈치 안전블록 · 구조장비

산업안전 보호구

전신형 안전대 · 랜야드 · 안전블록
작업 환경에 맞춘 개인보호구

산업용 소재·건축자재

산업용 테이프 · 접착제 · 기능성 소재
구조용 스크류 · 목조건축 자재

SNAPSHOOT
ENGINEERED FOR A SAFER TOMORROW

홈페이지 바로가기 →
문서 라이브러리 바로가기 →
시공사례 바로가기 →
문의하기 →

본 페이지는 스냅슛의 자료 제공 안내용 표지입니다. 첨부된 인증서 및 기술자료의 발행 주체, 인증 대상과 적용 범위는 해당 원문을 기준으로 확인해 주세요.

NEOGEN

Product Instructions
SESUP001

EN *Salmonella* Enrichment Supplement

Issue Date: 2024-01

NEOGEN®

Product Instructions

Salmonella Enrichment Supplement

DESCRIPTION

For use with Neogen® *Salmonella* Enrichment Base in conjunction with the Neogen® Petrifilm® *Salmonella* Express (SALX) System.

SAFETY

The user should read, understand, and follow all safety information in the instructions for the Neogen *Salmonella* Enrichment Supplement. Retain the safety instructions for future reference.

⚠ **WARNING** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

NOTICE Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

⚠ WARNING

To reduce the risks associated with false-negative results leading to release of contaminated food or beverage product:

- Always follow the proper enrichment scheme (time, temperature, and dilution).
- Always use care when pipetting or transferring samples.
- Follow the protocol and perform the tests exactly as stated in the product instructions for use.

To reduce the risks associated with false negative or false positive results leading to the release of contaminated product or the destruction of good product:

- Always follow product instructions for preparation, handling, and storage.

To reduce the risks associated with exposure to chemicals and biohazards:

- Dispose of enriched samples according to all applicable government regulations and applicable laboratory procedures.
- Read, understand, and follow the safety information in the Neogen Petrifilm *Salmonella* Express System Instructions for Use and the associated MSDS.
- Always follow standard good laboratory safety practices (GLP)[1], including proper containment procedures, wearing appropriate protective apparel and eye protection while handling testing materials and test samples.

To reduce the risks associated with exposure to chemicals which, if not avoided, may result in allergic skin reaction, allergic respiratory reaction, target organ effects, genotoxic or mutagenic effects and may interact with genetic material and possibly alter gene expression:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using this product.
- If eye contact does occur, flush eyes with large amounts of water. If signs/symptoms persist, get medical attention.
- If skin contact occurs, remove contaminated clothing and shoes. Immediately flush skin with large amounts of water. Get medical attention, wash contaminated clothing, and clean shoes before reuse.
- Do not eat, drink, or smoke when using this product.
- Always wash exposed areas thoroughly with soap and water.
- If inhaled, remove person to fresh air. If signs/symptoms develop, get medical attention.
- If swallowed, do not induce vomiting unless instructed to do so by medical personnel. Give victim two glasses of water. Never give anything by mouth to an unconscious person. Get medical attention.

To reduce the risks associated with environmental contamination:

- Follow current industry standards and local regulations for disposal of contaminated waste.

NOTICE: To reduce false-positive results due to cross contamination that result in re-testing or the rejection of food or beverage product:

- Always wear gloves.
- Use a new pipette tip for each sample transfer.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at **www.neogen.com,** or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

LIMITATION OF NEOGEN LIABILITY

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

STORAGE AND DISPOSAL

Storage: 2-8°C. Material is hygroscopic. Keep container tightly closed. Protect from light, heat and moisture.

Disposal: After use, enriched samples may contain microorganisms that may be a potential biohazard.

Follow current industry standards and local regulations for disposal of contaminated waste. Consult the Material Safety Data Sheet for additional information.

INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Wear appropriate protective apparel and follow standard good laboratory safety practices (GLP)[1].

Use an analytical balance capable of reading at least three to four decimal places to weigh the supplement.

Preparation of Supplement: Allow vial to come to ambient temperature. Gently invert the vial 3-4 times to mix the supplement powder. Aseptically add 50 mg (± 2.5 mg) of the Neogen® *Salmonella* Enrichment Supplement to 1 L of sterilized Neogen® *Salmonella* Enrichment Base tempered to 20 to 42.5°C. Mix thoroughly.

Composition: (typical formula g/L):

Selective Agents 0.05 g

REFERENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.

EXPLANATION OF SYMBOLS

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation
620 Lesher Place
Lansing, MI 48912 USA
neogen.com

Neogen Europe Ltd.
The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.
Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland



Neogen Corporation
620 Lesher Place Lansing, MI 48912 USA
www.neogen.com

Neogen is a trademark of Neogen Corporation.
© Neogen Corporation 2024.
All rights reserved.
FS00832A