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NEOGEN®

Meet the New *Listeria* Right Now™

Your environmental testing just got an upgrade.

Now integrated with the Neogen® Molecular Detection System™, the AOAC® PTM℠-validated *Listeria* Right Now™ features improved performance and a more efficient protocol. Make faster critical decisions with real-time results from user-friendly software.

Learn more at: **info.neogen.com/ListeriaRightNow**

Robust Performance

- **Detection where it counts.** Collect samples on a range of surfaces, like sealed concrete, stainless steel, and plastic
- **Tolerant to a range of sanitizers.** Chlorine-based, peroxyacetic acid, high acid cleaners, and quaternary ammonium
- **Results in about 2 hours.** That includes samples to actionable results
- **Unique sampler.** Can sample up to a 12” x 12” surface area
- **Same trusted sensitivity.** Can detect 2-5 CFUs
- **AOAC PTM validation.** AOAC certified for confident results

Neogen® Molecular Detection System

- User-friendly software with intuitive user interface for quick and easy setup
- Real-time results help make critical decisions faster
- Color-coded assays by organism and matching software screens allow for better tracking and easy identification

Streamlined Workflow

ANSR	Molecular Detection System
Prepare lysis tubes *after sample collection*	**Prepare lysis tubes** *after sample collection*
Transfer to heat block *37°C for 10 minutes*	**Transfer to heat block** *100 °C for 10 minutes*
Transfer in heat block *80°C for 20 minutes*	
Preheat reagent tubes *56°C for 3-5 minutes*	**Rest sample tubes in chill block** *at room temperature for 5 minutes*
Pipette 50 µL *of sample*	**Pipette 50 µL** *of sample into reagent tubes*
Transfer pipetted sample *into reagent tubes*	**Reconstitute the reagent** *pipette up and down five times*
Vortex *reagent sample tubes*	
Tap reagent tubes *to remove bubbles*	**Load reagent tubes** *into instrument for results*
Load reagent tubes *into instrument for results*	
9 STEPS	**6 STEPS**

Ordering Information

Description	Cat. No.	SKU	Contents
Neogen® Molecular Detection Assay – *Listeria* Right Now™	MDALRN96	700007660	96 Tests/Kit
Neogen® Molecular Detection Right Now™ Matrix Control	RNMC96	700007661	96 Tests/Kit
Neogen® Molecular Detection Instrument and Accessory Kit	MDS100	700002195	1 Each
Neogen® Molecular Detection Heat/Chill Block Insert	MDSHBIN	700000145	1 Each

Request a free trial at: **info.neogen.com/ListeriaRightNow**

NEOGEN®

PERFORMANCE TESTED AOAC RESEARCH INSTITUTE LICENSE NUMBER 042604

Neogen Corporation, 620 Lesher Place, Lansing, MI 48912 USA.
© Neogen Corporation 2026. All rights reserved. Neogen is a registered trademark of Neogen Corporation.
All other marks mentioned are the property of their respective owners.
FS02247_0326